

# ET 53-85

*Elevating bucky table*

CE



## Main features

- Floating patient table top with adjustable height
- Foot pedals control for up/down motorized motion and table top electromagnetic brakes
- Undertable bucky carriage movable by hand grip with push button control of electromagnetic brake
- Max. admitted weight for the patient of 150 kg
- Remote control of the table top up/down motion

## Features

<b>Bucky table</b>	2qws
• Patient table top	plastic laminated sheet
• Absorption according EN 60 601-2-35	0,9 mm Al eq.
• Longitudinal travel (without anticollision device)	110 cm (55 cm each side)
• Longitudinal travel (with anticollision device)	90 cm (45 cm each side)
• Transversal travel	26 cm (+/- 13 cm)
• Control of electromagnetic brakes	by foot pedals
• Table top to film distance	7 cm
• Table top to floor minimum distance	53 cm
• Table top height variation	32 cm
• Average speed	3 cm/s
• Patient max weight allowed	150 kg
<b>Bucky carriage</b>	
• Longitudinal travel	50,5 cm
• Electromagnetic brakes control	by push button
<b>Power requirements</b>	
• Line voltage / frequency	230 (220/240) Vac 50/60 Hz
• Max. absorbed power	300 VA
<b>Standard color</b>	
• Painted parts	RAL 9002
• Adhesive overlays	Gray
<b>Weights</b>	
• Elevating bucky table	280 kg
<b>Equipment identification</b>	
• Classification	class I - type B
• Tipology of installation	permanently connected
• Working mode	continous with intermittent load
<b>Standard compliance</b>	
• Safety	EN 60 601
• Electromagnetic compatibility	89/336 European directive
<b>Transport and stock environment</b>	
• Temperature range	- 15 °C ÷ 60 °C
• Relative humidity range (without condensing)	10 % ÷ 90 %
• Atmospheric pressure range	500 hPa ÷ 1060 hPa
<b>Operating environment</b>	
• Temperature range	10 °C ÷ 40 °C
• Relative humidity range (without condensing)	30 % ÷ 75 %
• Atmospheric pressure range	700 hPa ÷ 1060 hPa

### Accessories to complete the table.


PB0LRE	Bucky / cassette tray
KGB104	Grid K FFD =100 cm R=10:1 40 l/cm

### Accessories on request

ZA0012	Pair of patient hand grips	ZA0013	24x30 cm articulated lateral cassette holder
ZA0011	Compression band with ratchet	ZA0014	35x43 cm lateral cassette holder
ZA0015	Pair of leg support	CTA004	Flebo support
ETA002	Table top up/down remote control foot pedal	CTA020	Fixing clamp for flebo support

### Options

ETA001	Low absorption table top (0,45 mm Al eq.)
ETA003	Table top anti-collision device
	Table top positive brakes
	Bucky with automatic cassette loader
	A.E.C - Automatic Exposure Control system



Bucky table with height adjustment cod. ET8420